



YIPF5 Antibody

Product Code	CSB-PA026261GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969M3
Immunogen	Human YIPF5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Yip1 domain family, member 5;YIPF5;DKFZp313L2216;FinGER5;SB140;SMAP-5;SMAP5;YIP1A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	YIPF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.