



# YME1L1 Antibody

<b>Product Code</b>	CSB-PA026267GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96TA2
<b>Immunogen</b>	Human YME1L1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	ELISA, WB, IHC
<b>Form</b>	Liquid
<b>Storage Buffer</b>	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	YME1L1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.