



ZBTB7B Antibody

Product Code	CSB-PA026344GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15156
Immunogen	Human ZBTB7B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	zinc finger and BTB domain containing 7B;ZBTB7B;DKFZp686G01254;THPOK;ZBTB15;ZFP67;ZNF857B;c-Krox;hcKrox ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZBTB7B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.