



# ZC3HC1 Antibody

|                            |                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA026371GA01HU                                                                   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                        |
| <b>Uniprot No.</b>         | Q86WB0                                                                               |
| <b>Immunogen</b>           | Human ZC3HC1                                                                         |
| <b>Raised In</b>           | Rabbit                                                                               |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                      |
| <b>Tested Applications</b> | ELISA,WB                                                                             |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Purification Method</b> | Antigen Affinity Purified                                                            |
| <b>Isotype</b>             | IgG                                                                                  |
| <b>Alias</b>               | zinc finger, C3HC-type containing 1;ZC3HC1;1110054L24Rik;HSPC216;NIPA ;              |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                       |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                 |
| <b>Target Names</b>        | ZC3HC1                                                                               |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.          |