



ZDHHC5 Antibody

Product Code	CSB-PA026419GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0B5
Immunogen	Human ZDHHC5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	zinc finger, DHHC-type containing 5;ZDHHC5;DKFZp586K0524;KIAA1748;ZNF375 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZDHHC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.