



ZBTB18 Antibody

Product Code	CSB-PA026604GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99592
Immunogen	Human ZNF238
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	zinc finger protein 238;ZNF238;C2H2-171;RP58;TAZ-1;ZBTB18 ;ZFP238
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZBTB18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.