



MTDH Antibody

Product Code	CSB-PA015143GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UE4
Immunogen	Human AEG-1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	metadherin;MTDH;3D3;AEG-1;AEG1;LYRIC;LYRIC/3D3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MTDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.