



ARHGDI Antibody

Product Code	CSB-PA002037GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52565
Immunogen	Human ARHGDI,aGDI
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Protein A purified
Isotype	IgG
Alias	Rho GDP dissociation inhibitor (GDI) alpha;ARHGDI;GDIA1;MGC117248;RHOGDI;RHOGDI-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGDI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.