



CTNNB1 Antibody

Product Code	CSB-PA006169GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Human b-cat
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	catenin (cadherin-associated protein), beta 1, 88kDa;CTNNB1;CTNNB;DKFZp686D02253;FLJ25606;FLJ37923 ;?? catenin;BETA catenin;CTNNB1
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Adhesion
Target Names	CTNNB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.