



# Recombinant Human Serine/threonine-protein phosphatase 5 (PPP5C)

<b>Product Code</b>	CSB-BP018582HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53041
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AMAEGERTE CAEPPRDEPP ADGALKRAEE LKTQANDYFK AKDYENAIKF YSQAIELNPS NAIYYGNRSL AYLRTCEYGY ALGDATRAIE LDKKYIKGYY RRAASNALG KFRAALRDYE TVVKVPHDK DAKMKYQECN KIVKQKA FER AIAGDEHKRS VVDSLDIISM TIEDEYSGPK LEDGKVTISF MKELMQWYKD QKKLHRKCA Y QILVQVKEVL SKLSTLVETT LKETEKITVC GDTHGQFYDL LNIFELNGLP SETNPYIFNG DFVDRGSFSV EVILTLFGFK LLYPDHFHLL RGNHETDNMN QIYGFEDEVK AKYTAQMYEL FSEVFEWLPL AQCINGKVL I MHGGLFSEDG VTLDDIRKIE RNRQPPDSGP MCDLLWSDPQ PQNGRSISK R GVSCQFGPDV TKAFLEENL DYIIRSHEVK AEGYEVAHGG RCVTVFSAPN YCDQMGNKAS YIHLQGSCLR PQFHQFTAVP HPNVKPMAYA NTLQLGMM
<b>Source</b>	Baculovirus
<b>Target Names</b>	PPP5C
<b>Protein Names</b>	Recommended name: Serine/threonine-protein phosphatase 5 Short name= PP5 EC= 3.1.3.16 Alternative name(s): Protein phosphatase T Short name= PP-T Short name= PPT
<b>Expression Region</b>	2-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.