



NCR2 Antibody

Product Code	CSB-PA015550GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95944
Immunogen	Human CD336
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	natural cytotoxicity triggering receptor 2;NCR2;CD336;LY95;NK-p44;NKP44;dJ149M18.1 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.