



Recombinant Danio rerio Coiled-coil domain-containing protein 69 (ccdc69)

Product Code	CSB-BP004731DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7MC22
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	GCHNSKVCG QVSKKKKKNK AQEGEKRTKD IKLQDGRGTH SSEEERCCLN NQLENYEWQL KILHAVLTAS GDQEREQLLK DHPGDICTLV HSITEKVKTE ISADLNDLHE QQMRSVSEQH QSETEELQRL HSEEKNVLE SHAAEEALK DQIEDLTSEL KLFNELKRRR EASTLKRDLQ RNIETHGSPG EFWEQEQESL LFVIEMKRQN LQDQGNKLLQ MEAQVEKNLS LEDQLLQALQ QNEDFRVRIE NYQSLIQQLS KEQNDMQEVL EKQSLQNQKL IQEKEELLFK LLHRRDSCST FHLPSVIPTQ VSPS
Source	Baculovirus
Target Names	ccdc69
Protein Names	Recommended name: Coiled-coil domain-containing protein 69
Expression Region	2-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.