



MARK3 Antibody

Product Code	CSB-PA013498GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27448
Immunogen	Human CTAK1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	MAP/microtubule affinity-regulating kinase 3;MARK3;CTAK1;KP78;PAR1A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MARK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.