



IL13RA2 Antibody

Product Code	CSB-PA011592GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14627
Immunogen	Human IL13RA2,CD213A2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	interleukin 13 receptor, alpha 2;IL13RA2;CD213A2;CT19;IL-13R;IL13BP ;IL13Ralpha2; IL-13R-alpha-2
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL13RA2
Target Details	This protein is closely related to IL13RA1, a subunit of the interleukin 13 receptor complex. This protein binds IL13 with high affinity, but lacks cytoplasmic domain, and does not appear to function as a signal mediator. It is reported to play a role in the internalization of IL13.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.