



MAS1L Antibody

Product Code	CSB-PA013506GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35410
Immunogen	Human MAS1L-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	MAS1 oncogene-like;MAS1L;MAS-L;MGC119987;MRG;dJ994E9.2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAS1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.