



MEF2B Antibody

Product Code	CSB-PA013672GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02080
Immunogen	Human MEF2B-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	LOC729991-MEF2B readthrough transcript; LOC729991-MEF2B; MEF2B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MEF2B
Target Details	The product of this gene is a member of the MADS/MEF2 family of DNA binding proteins. The protein is thought to regulate gene expression, including expression of the smooth muscle myosin heavy chain gene. This region undergoes considerable alternative splicing, with transcripts supporting two non-overlapping loci (geneID:729991 and 100271849) as well as numerous read-through transcripts that span both loci (annotated as geneID:4207). Several isoforms of this protein are expressed from either this locus or from some of the read-through transcripts annotated on geneID:4207.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.