



CISD1 Antibody

Product Code	CSB-PA005442GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ45
Immunogen	Human mitoNEET,CISD1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	CDGSH iron sulfur domain 1;CISD1;C10orf70;MDS029;MGC14684;ZCD1;mitoNEET ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CISD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.