



N6AMT1 Antibody

Product Code	CSB-PA010286GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5N5
Immunogen	Human N6AMT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	N-6 adenine-specific DNA methyltransferase 1 (putative);N6AMT1;C21orf127;HEMK2;MGC19995;MTQ2;N6AMT;PRED28 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	N6AMT1
Target Details	This protein belongs to the methyltransferase superfamily. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.