



POU5F1 Antibody

Product Code	CSB-PA018403GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01860
Immunogen	Human OCT4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	POU class 5 homeobox 1;POU5F1;MGC22487;OCT3;OCT4;OTF3;OTF4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POU5F1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.