



PACS2 Antibody

Product Code	CSB-PA017372GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VP3
Immunogen	Human PACS2-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phosphofurin acidic cluster sorting protein 2;PACS2;FLJ25488;KIAA0602;PACS1L ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PACS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.