



PPP1R11 Antibody

Product Code	CSB-PA018508GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60927
Immunogen	Human PPP1R11-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	protein phosphatase 1, regulatory (inhibitor) subunit 11;PPP1R11;HCG-V;HCGV;MGC125741;MGC125742;MGC125743;TCTE5;TCTEX5 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PPP1R11
Target Details	This gene encodes a specific inhibitor of protein phosphatase-1 (PP1) with a differential sensitivity toward the metal-independent and metal-dependent forms of PP1. The gene is located within the major histocompatibility complex class I region on chromosome 6.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.