



IL18 Antibody

Product Code	CSB-PA011608GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14116
Immunogen	Human Pre-IL18
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC, IF
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	interleukin 18 (interferon-gamma-inducing factor); IL18; IGIF; IL-18; IL-1g; IL1F4; MGC12320 ;
Immunogen Species	Homo sapiens (Human)
Target Names	IL18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.