



ATP2A2 Antibody

Product Code	CSB-PA002333GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16615
Immunogen	Human SERCA2,ATP2A2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP2A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.