



SLC2A1 Antibody

Product Code	CSB-PA021546GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11166
Immunogen	Human SLC2A1, GLUT1
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 2 (facilitated glucose transporter), member 1; SLC2A1; DYT17; DYT18; GLUT; GLUT1; MGC141895; MGC141896; PED ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC2A1
Target Details	This gene encodes a major glucose transporter in the mammalian blood-brain barrier. Mutations in this gene have been found in a family with paroxysmal exertion-induced dyskinesia.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.