



SMYD1 Antibody

Product Code	CSB-PA021862GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NB12
Immunogen	Human SMYD1-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SET and MYND domain containing 1;SMYD1;BOP;ZMYND18;ZMYND22 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMYD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.