



SRXN1 Antibody

Product Code	CSB-PA022698GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYN0
Immunogen	Human SRX1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sulfiredoxin 1 homolog (S. cerevisiae);SRXN1;C20orf139;Npn3;SRX1;YKL086W;dJ850E9.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SRXN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.