



TFAP2A Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA023414GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P05549 |
| Immunogen | Human TFAP2A,AP-2 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | transcription factor AP-2 alpha (activating enhancer binding protein 2 alpha);TFAP2A;AP-2;AP-2alpha;AP2TF;BOFS;FLJ51761;TFAP2 ;AP2 |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TFAP2A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |