



TINF2 Antibody

Product Code	CSB-PA023566GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSI4
Immunogen	Human TIN2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	TERF1 (TRF1)-interacting nuclear factor 2;TINF2;TIN2;TIN2L ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TINF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.