



UNC5D Antibody

Product Code	CSB-PA025632GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXZ4
Immunogen	Human UNC5D-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	unc-5 homolog D (C. elegans);UNC5D;FLJ16019;KIAA1777;PRO34692;Unc5h4 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UNC5D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.