



SLC18A2 Antibody

Product Code	CSB-PA021427GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05940
Immunogen	Human VMAT2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 18 (vesicular monoamine), member 2; SLC18A2; MGC120477; MGC120478; MGC26538; SVAT; SVMT; VAT2; VMAT2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC18A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.