



# Recombinant *Shewanella pealeana* Queuine tRNA-ribosyltransferase (tgt)

<b>Product Code</b>	CSB-YP019146STP
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A8H2K8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Shewanella pealeana</i> (strain ATCC 700345 / ANG-SQ1)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKFELDTTQG RARRGRLIFE RGTVETPAFM PVGTYGTVKG MTPPEEVRETG ADILLGNTFH LWLRPGEEIM RKHGDLHDFM NWQRPILTDS GGFQVFLSGD IRKITEEGVH FRSPINGEKI FLDPEKSMQI QDSLGS DVVM IFDECTPYPA TEDEARKSMQ MSLRWAQRSR DEFDRLENPN SLFGIIQGGV YEDLRDESLN GLVDIGFDGY AVGGLAVGEP KEDMHRILEH VCPKIPADKP RYLMGVGKPE DLVEGVRRGV DMFDCVMPTR NARNGHLFTS EGVKIRNAR HRDDTSPLDD KDCYTCCKNY SRAYLYHLDR CNEILGARLN TIHNLRY YQR LMEGLRGAIE TGTLDAFVTE FYTSQGREVP EVPELSD
<b>Source</b>	Yeast
<b>Target Names</b>	tgt
<b>Protein Names</b>	Recommended name: Queuine tRNA-ribosyltransferase EC= 2.4.2.29 Alternative name(s): Guanine insertion enzyme tRNA-guanine transglycosylase
<b>Expression Region</b>	1-377
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.