



Recombinant Human Tumor protein D52 (TPD52)

Product Code	CSB-BP024096HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55327
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDCREMDLYE DYQSPDFDA GVNKSYLYLS PSGNSSPPGS PTLQKFGLLR TDPVPEEGED VAATISATET LSEEEQEELR RELAKVEEEI QTLSQVLA EKHLAEIKRK LGINSLQELK QNIAKGWQDV TATSAYKKT ETLSQAGQKA SAAFSSVGSV ITKKLEDVKN SPTFKSFEEK VENLKSXVGG TKPAGGDFGE VLNSAANASA TTTEPLPEKT QESL
Source	Baculovirus
Target Names	TPD52
Protein Names	Recommended name: Tumor protein D52 Alternative name(s): Protein N8
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.