



# Recombinant Danio rerio RUN domain-containing protein 3A (rundc3a)

<b>Product Code</b>	CSB-BP020590DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A7YDW0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEQSFMQSAM AMGAQSKKGF SRSIAVERKN LITVCRFSVK TLLEKYTAEP IDDSSEEFIN FAAILEHILS HRFKGGSGWF DGQRSFWDPI RLACSKVPNN CISSIENMEN INSSRAKGRA WLRVALMEKR LSEYIATALR DSRTTRRFYD EGAIMLREEA TVLTGMLIGL GAIDFSFCLK GEALDGKSEA VIDYTPYLKF TQSYDYLSDD EDGQSVSSN SDDSVHEPYI PLVTDEESWR NKCRKMEQRF KIVNAQKGYL EELVRLRESQ LKNVEMENKK LTAGLEELQQ QSQKEKRELE NIILELQEQL TSLIPGESHP LSKDLSIPLV NQWPSLQNYN NQEDSKLYHR GSFSPPEPHI SLTTGSQRTE RKQNGKSWCT AEKDYTPSML GLCGSMSSLP SCKSLPSLRS TECLVNISAE PSPALTPS
<b>Source</b>	Baculovirus
<b>Target Names</b>	rundc3a
<b>Protein Names</b>	Recommended name: RUN domain-containing protein 3A Alternative name(s): Rap2-interacting protein Short name= Rap2ip
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.