



Recombinant Danio rerio RUN domain-containing protein 3A (rundc3a)

Product Code	CSB-YP020590DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7YDW0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEQSFMQSAM AMGAQSKKGF SRSIAVERKN LITVCRFSVK TLLEKYTAEP IDDSSEEFIN FAAILEHILS HRFKGGSGWF DGQRSFWDPI RLACSKVPNN CISSIENMEN INSSRAKGRA WLRVALMEKR LSEYIATALR DSRTTRRFYD EGAIMLREEA TVLTGMLIGL GAIDFSFCLK GEALDGKSEA VIDYTPYLKF TQSYDYLSDD EDGQSVSSN SDDSVHEPYI PLVTDEESWR NKCRKMEQRF KIVNAQKGYL EELVRLRESQ LKNVEMENKK LTAGLEELQQ QSQKEKRELE NIILELQEQL TSLIPGESHP LSKDLSIPLV NQWPSLQNYN NQEDSKLYHR GSFSPPEPHI SLTTGSQRTE RKQNGKSWCT AEKDYTPSML GLCGSMSSLP SCKSLPSLRS TECLVNISAE PSPALTPS
Source	Yeast
Target Names	rundc3a
Protein Names	Recommended name: RUN domain-containing protein 3A Alternative name(s): Rap2-interacting protein Short name= Rap2ip
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.