



# Recombinant Human Zinc finger protein 131 (ZNF131)

<b>Product Code</b>	CSB-BP026530HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P52739
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEAEETMECLQEFPEHHKMILDRLENEQREQDRFTDITLIVDGHFFKAHKAVLA ACSKFFY KFFQEFTQEPLVEIEGVSKMAFRHLIEFTYTAKLMIQGEEEEANDVWKAAEFLQM LEAIKA LEVRNKENSAPLEENTTGKNEAKKRKIAETSNVITESLPSAESEPVEIEVEIAEG TIEVE DEGIETLEEVASAKQSVKYIQSTGSSDDALALLADITSKYRQGDRKGQIKEDG CPSDPT SKQVEGIEIVELQLSHVKDLFHCEKCNRSFKLFYHFKEHMKSHSTESFKCEICN KRYLRE SAWKQHLNLCYHLEEGGVSKKQRTGKKIHVCQYCEKQFDHFGHFKEHLRKHT GEKPFECPN CHERFARNSTLKCHLTACQTGVGAKKGRKKLYECQVCNSVFNSWDQFKDHL VIHTGDKPN HCTLCDLWFMQGNELRRHLSDAHNIERLVTEEVLSVETRVQTEPVTSMTIIEQ VGKVHV LPLLQVQVDSAQVTVEQVHPDLLQDSQVHDSHMSELPEQVQVSYLEVGRIQT EEGTEVHV EELHVERVNQMPVEVQTELLEADLDHVTPEIMNQEERESSQADAAEAAREDH EDAEDLET KPTVDSEAEKAENEDRTALPVLE
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF131
<b>Protein Names</b>	Recommended name: Zinc finger protein 131
<b>Expression Region</b>	1-623
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.