



Recombinant Human Zinc finger protein 131 (ZNF131)

Product Code	CSB-EP026530HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P52739
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MEAEETMECLQEFPEHHKMILDRLENEQREQDRFTDITLIVDGHFFKAHKAVLA ACSKFFY KFFQEFTQEPLVEIEGVSKMAFRHLIEFTYTAKLMIQGEEEEANDVWKAAEFLQM LEAIKA LEVRNKENSAPLEENTTGKNEAKKRKIAETSNVITESLPSAESEPVEIEVEIAEG TIEVE DEGIETLEEVASAKQSVKYIQSTGSSDDALALLADITSKYRQGDRKGQIKEDG CPSDPT SKQVEGIEIVELQLSHVKDLFHCEKCNRSFKLFYHFKEHMKSHSTESFKCEICN KRYLRE SAWKQHLNLCYHLEEGGVSKKQRTGKKIHVCQYCEKQFDHFGHFKEHLRKHT GEKPFECPN CHERFARNSTLKCHLTACQTGVGAKKGRKKLYECQVCNSVFNSWDQFKDHL VIHTGDKPN HCTLCDLWFMQGNELRRHLSDAHNISRVLVTEEVLSVETRVQTEPVTSMTIIEQ VGKVHV LPLLQVQVDSAQVTVEQVHPDLLQDSQVHDSHMSELPEQVQVSYLEVGRIQT EEGTEVHV EELHVERVNQMPVEVQTELLEADLDHVTPEIMNQEERESSQADAAEAAREDH EDAEDLET KPTVDSEAEKAENEDRTALPVLE</p>
Source	E.coli
Target Names	ZNF131
Protein Names	Recommended name: Zinc finger protein 131
Expression Region	1-623
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.