



Recombinant Danio rerio Lysine-specific demethylase 8 (jmjd5)

Product Code	CSB-MP011953DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8E534
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MASVWTDIRA VLPSTVSEFP LDFSEKIDLS VLKCLELSRD QLYSEADCPV SAERAQIIID YSWEKLNIGT WRDVDKEWRR VYSYGCLFKV LSLCHGNPPQ NIIQEAVRTC DMSLLMGAAI MDNILQRLVG ILRNKIKTTC PNKAERSEEP FSKKRKHDK SEPVLNPTKE VPRIHCPSLE RFRSDFLDPK KPVIIEGID LWPAFTQHPW SIDYLRTVAG CRTVPIEVGS KYTDEEWSQK LITVNDFIDR YITGTEEDGV GYLAQHQLFD QVPELKEDIR IPDYCCLGEG DEDDITINAW FGPGGTVSPL HQDPQQNFLA QVVGRKYIRL YSPEDTKSLY PHESQLLHNT SQVEVENPDL VKFPDFSRAS YEECVLCPGD VLFIPLQHWY YVRSLELSFS VSFWWS
Source	Mammalian cell
Target Names	kdm8
Protein Names	Recommended name: Lysine-specific demethylase 8 EC= 1.14.11.27 Alternative name(s): JmjC domain-containing protein 5 Jumonji domain-containing protein 5
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.