



Recombinant Azorhizobium caulinodans Lipoyl synthase (lipA)

Product Code	CSB-BP012927DPR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8I4L8
Product Type	Recombinant Protein
Immunogen Species	Azorhizobium caulinodans (strain ATCC 43989 / DSM 5975 / JCM 20966 / NBRC 14845 / NCIMB 13405 / ORS 571)
Purity	>85% (SDS-PAGE)
Sequence	MVTVVNTLNR PRHPEKQNRP ETEVLRKPDW IRVKAPGSAG WSNTAGIVRA NGLHTVCEEA GCPNIGECWE KKHATFMIMG DTCTRACSFC NVRTGMPKAL DLDEPQKVGE AVAKLGLSHV VITSVDRDDL TDGGAEHFAR TIASIRKLSP GTTIEILTPD FLRKPGEAIEV VVAARPDVFN HNLETVPGKY LTVRPGARYF HSLRLLQQVK ELDPSIFTKS GIMVGLGEER NEVLQLMDDL RAAEVDFMTI GQYLQPTRKH HKVERFVTPD EFKAYETVAY AKGFLMVSSS PLTRSSHAG DDFAKLRAAR EAKLGRI
Source	Baculovirus
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.