



# Recombinant Chicken Replication factor C subunit 2 (RFC2)

<b>Product Code</b>	CSB-EP019589CH-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53033
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEEEEVLEVV EDEKAGPAAA EKRGPDTLGL SAPAASGHYE LPWVEKYRPL KLCEVVGNE DTVSRLEVF AKEGNVPNIIA GPPGTGKTTS ILCLARALLG PALKDAVLEL NASNDRGIDV VRNKIKMFAQ QKVTLPKGRH KIIILDEADS MTDGAQQALR RTMEIYSKTT RFALACNASD KIIIEPIQSRC AVLRYTKLTD SQILARLLKI VEKEDVPYTD DGLEAIIFTA QGDMRQALNN LQSTYSGFGF INSENVFKVC DEHPPLLVE MIQHCINANI DEAYKIL A HLWRLGYSPEDV IGNIFRVCKT FQMP EY LKLE FIK EIGYTHM KIAEGVNSLL QMAGLLARLC QKTAAPAAS
<b>Source</b>	E.coli
<b>Target Names</b>	RFC2
<b>Protein Names</b>	Recommended name: Replication factor C subunit 2 Alternative name(s): Activator 1 40 kDa subunit Short name= A1 40 kDa subunit Activator 1 subunit 2 Replication factor C 40 kDa subunit Short name= RF-C 40 kDa subunit Shor
<b>Expression Region</b>	1-359
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the action of the accessory proteins, proliferating cell nuclear antigen (PCNA) and replication factor C (RFC). RFC, also called activator 1, is a protein complex consisting of five distinct subunits of 145, 40, 38, 37, and 36.5 kD. This gene encodes the 40 kD subunit, which has been shown to be responsible for binding ATP. Deletion of this gene has been associated with Williams syndrome. Alternatively spliced transcript variants encoding distinct isoforms have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.