



Recombinant Escherichia coli O9:H4 Phosphoserine aminotransferase (serC)

Product Code	CSB-EP018838EJF-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZYL0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAQIFNFSSG PAMPLAEVLK QAQQELRDWN GLGTSVMEVS HRGKEFIQVA EEAEKDFRDL LNVPSNYKVL FCHGGGRGQF AAVPLNILGD KTTADYVDAG YWAASAIKEA KKYCTPNVFD AKVTVDGLRA VKPMREWQLS DNAAYMHYCP NETIDGIAID ETPDFGKDVV VAADFSSTIL SRPIDVSRYG VIYAGAQKNI GPAGLTIVIV REDLLGKANI ACPSILDYSI LNDNGSMFNT PPTFAWYLSG LVFKWLKANG GVAEMDKINQ QKAELLYGVI DNSDFYRNDV AKANRSRMNV PFQLADSALD KLFLEESFAA GLHALKGHRV VGGMRASIYN AMPLEGVKAL TDFMVEFERR HG</p>
Source	E.coli
Target Names	serC
Protein Names	<p>Recommended name: Phosphoserine aminotransferase EC= 2.6.1.52 Alternative name(s): Phosphohydroxythreonine aminotransferase Short name= PSAT</p>
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>