



Recombinant *Drosophila melanogaster* Aprataxin and PNK-like factor (CG6171)

Product Code	CSB-BP001908DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8JR14
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSATDASTAD SGAKRKSSSED ITHNCNANFG AENGLRKRVK SEEPVASIKD ETNPEVPMKI KAEPVENADE PTSTTPAIKI KAEPADNGNS PAAAMVKTEP TNSNAQDAAD ESTVSSSSIR TSCRFGIRCY RRNPAHRSAE AHPGDQDYRR PNFPAPPLGT PACPFGNACY RRNPVHFQDY SHPADFNSAQ NIRNRLRQRR AQRQNDDDDSG TDEEDEPFGG DNDRDADYRP GADINEDEDD ELEFDSQPIS GDDYD
Source	Baculovirus
Target Names	CG6171
Protein Names	Recommended name: Aprataxin and PNK-like factor EC= 4.2.99.18 Alternative name(s): Apurinic-aprimidinic endonuclease APLF
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.