



# Recombinant Mouse Voltage-dependent L-type calcium channel subunit beta-3 (Cacnb3)

<b>Product Code</b>	CSB-YP004413MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P54285
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MYDDSYVPGF EDSEAGSADS YTSRPSLDSD VSLEEDRESA RREVESQAQQ  QLERAKHKPV AFAVRTNVS Y CGVLDEECPV QGSGVNFEAK DFLHIKEKYS  NDWWIGRLVK EGGDIAFIPS PQRLESIRLK QEQAARRSGN PSSLGDIGNR  RSPPPSLAKQ KQKQAEHVPP YDVVPSMRPV VLVGPSLKGY EVTDMMQKAL  FDFLKHFRFDG RISITRV TAD LSLAKRSVLN NPGKRTIER SSARSSIAEV  QSEIERIFEL AKSLQLVLD ADTINHPAQL AKTSLAPIIV FVKVSSPKVL  QRLIRSRGKS QMKHLTVQMM AYDKLVQCPP ESFDVILDEN QLEDACEHLA  EYLEVYWRAT HHPAPGPGLL GPPSAIPGLQ NQQLLGERVE EHSPLERDSL  MPSDEASESS RQAWTGSSQR SSRHLEEDYA DAYQDLYQPH  RQHTSGLPSA NGHDPQDRLL AQDSEHDHND RNWQRNRPWP KDSY</p>
<b>Source</b>	Yeast
<b>Target Names</b>	Cacnb3
<b>Protein Names</b>	Recommended name: Voltage-dependent L-type calcium channel subunit beta-3 Short name= CAB3 Alternative name(s): Calcium channel voltage-dependent subunit beta 3 Short name= CCHB3
<b>Expression Region</b>	1-484
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.