



# Recombinant Rat Mothers against decapentaplegic homolog 3 (Smad3)

<b>Product Code</b>	CSB-MP021788RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P84025
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>SSILPFTPP IVKRLLGWKK GEQNGQEEKW CEKAVKSLVK KLKKTGQLDE            LEKAITQNV NTKCITIPRS LDGRLQVSHR KGLPHVIYCR LWRWPDLHSH            HELRAMELCE FAFNMKKDEV CVNPYHYQRV ETPVLPPVLV PRHTEIPAEF            PPLDDYSHSI PENTNFPAGI EPQSNIPETP PPGYLSEGE TSDHQMNHSM            DAGSPNLSPN PMSPAHHNLD LQPVTYCEPA FWCSISYYEL NQRVGETFHA            SQPSMTVDGF TDPSNSERFC LGLLSNVNRN AAVELTRRHI GRGVRLYYIG            GEVFAECLSD SAIFVQSPNC NQRYGWHPAT VCKIPPGCNL KIFNNQEFAA            LLAQSVNQQG EAVYQLTRMC TIRMSFVKGW GAEYRRQTVT STPCWIELHL            NGPLQWLDKV LTQMGSPSIR CSSVS</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Smad3
<b>Protein Names</b>	Recommended name: Mothers against decapentaplegic homolog 3 Short name= MAD homolog 3 Short name= Mad3 Short name= Mothers against DPP homolog 3 Alternative name(s): SMAD family member 3 Short name= SMAD 3 Short nam
<b>Expression Region</b>	2-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene mothers against decapentaplegic (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein functions as a transcriptional modulator activated by transforming growth factor-beta and is thought to play a role in the regulation of carcinogenesis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.