



# Recombinant Danio rerio Coiled-coil domain-containing protein 89 (ccdc89)

<b>Product Code</b>	CSB-BP004758DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A8WG43
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSPQKIPID LKMMIKDTKQ DMDDVHQALE KLRSLSQEER TEAEVLSRI DEQSTLICIL KRRADEMLLR CQALEKINAE LENLRANVQV ELENERKKSE LLEQRFMDLA SNHRELINFK NEYKSQNAKL EKENERLRKE NEALFSEEVQ EKEETILKLT QELKDLAEQH KILENEYQEK TTGFQIKLKE LMNLHRIKEA SLKNELSDTQ KQLKTAVEMC AELDLKLRQT QVKESMRETD TQERLEKLIK EKDELIDLSV QRGKIIRDKQ VEIQELERKR QEAENGRRHA EDRFEKEAAA VNSELRVKEL REALERA EFS CNELKKEFEA FKKHSCELLE KEKELNCKLR HMII
<b>Source</b>	Baculovirus
<b>Target Names</b>	ccdc89
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 89 Alternative name(s): Bc8 orange-interacting protein
<b>Expression Region</b>	1-354
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.