



RuBisCO activase Antibody, FITC conjugated

Product Code	CSB-PA00170G0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10871
Immunogen	Ribulose-1,5-bisphosphate carboxylase oxygenase (Native Protein)
Raised In	Rabbit
Species Reactivity	Spinach
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribulose bisphosphate carboxylase/oxygenase activase, chloroplastic (RA) (RuBisCO activase)
Immunogen Species	Spinach
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.