



TXN Antibody, Biotin conjugated

Product Code	CSB-PA02814D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10599
Immunogen	Recombinant Human Thioredoxin protein (2-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin (Trx) (ATL-derived factor) (ADF) (Surface-associated sulphydryl protein) (SASP), TXN, TRDX TRX TRX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TXN