



abscisic acid Monoclonal Antibody

Product Code	CSB-MA004211I0m
Storage	Shipped at 4°C Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	abscisic acid-BSA conjugate
Raised In	mouse
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG1
Clonality	monoclonal
Product Type	Small Molecules Antibody
Accession NO.	1A1D4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.