



# Ractopamine Monoclonal Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-MA001811I0m  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Immunogen</b>           | Ractopamine -BSA conjugate   |
| <b>Raised In</b>           | mouse  |
| <b>Specificity</b>         | No significant cross-reactivity or interference was observed                     |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Purification Method</b> | >95% protein G purified  |
| <b>Isotype</b>             | IgG1   |
| <b>Clonality</b>           | monoclonal   |
| <b>Product Type</b>        | Small Molecules Antibodies   |
| <b>Target Names</b>        | RA   |
| <b>Accession NO.</b>       | 5F3C2  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |