





HLA-DRB1 Antibody, Biotin conjugated

Product Code	CSB-PA17949D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01912
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DRB1-3 chain protein (31-266AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DRB1-3 chain (Clone P2-beta-3) (MHC class II antigen DRB1*3), HLA-DRB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HLA-DRB1

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA17949D0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA17949D0Rb at dilution of 1:100