



TGFB2 Antibody, Biotin conjugated

Product Code	CSB-PA07315D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61812
Immunogen	Recombinant Human Transforming growth factor beta-2 proprotein protein (304-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-2 (TGF-beta-2) (BSC-1 cell growth inhibitor) (Cetermin) (Glioblastoma-derived T-cell suppressor factor) (G-TSF) (Polyergin) [Cleaved into: Latency-associated peptide (LAP)], TGFB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TGFB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.